

Multicancer Early Detection (MCED)


Eric A. Klein, MD

Distinguished Scientist, GRAIL, Inc

Emeritus Professor and Chair, Glickman Urological and Kidney Institute
Cleveland Clinic Lerner College of Medicine

Globally, Cancer Cases Are On the Rise

Annual ASIRs for all cancers combined increased between 2007 and 2017 in 123 countries¹



ASIR, age-standardized incidence rate.

¹Fitzmaurice C et al. *JAMA Oncol.* 2019;5:1749-1768. ²Sung H et al. *CA Cancer J Clin.* 2021;71:209-249.

The Status Quo

Despite this:

USPSTF Recommendations for Cancer Screening

Cancer	Grade	Population	Modality/ Recommendation
Cervical	A	Women aged 21 to 65	Regular screening (3–5 years) using cervical cytology and/or HPV tests
Colorectal	A B	Adults aged 50 to 75 Adults aged 45–49	Regular annual screening, multiple effective methods available
Breast	B C	Women aged 50 to 74 Women aged 40 to 49	Biennial screening mammography
Lung	B	Adults aged 55–80, with history of smoking	Annual low-dose computed tomography (LDCT) screening
Prostate	C	Men aged 55 to 69	Periodic PSA screening on case-by-case basis




Mortality




> 600,000 people die of cancer every year in the US

Why is this Necessary?

Annual invasive cancer incidence in persons aged 50-79





Screened Patients are More Likely to be Diagnosed with Non-Screened Cancers




Paradigm Shift

Screening for individual cancers



Screening individuals for cancer(s)




- Breast cancer
- Lung cancer
- Colon cancer
- Prostate cancer
- Cervical cancer
- Lymphoid neoplasm
- Plasma-cell neoplasm
- Ovarian cancer
- Bladder cancer
- Gastrointestinal cancer
- Liver cancer
- Pancreatic cancer
- Head-and-neck cancer
- Anorectal cancer
- Uterine cancer
- Kidney cancer
- Melanoma
- Thyroid
- Myeloid neoplasm
- Sarcoma
- Multiple other cancers

GRAIL Clinical Development Program

Test Development, Validation, and Implementation in Population-Scale Studies

1	CCGA (n=15,254)	Develop and validate a cell-free DNA-based MCEd test <i>Enrollment: complete, published</i>	<i>Annals of Oncology and Cancer Cell 2020-2023</i>
2	PATHFINDER (n=6,662)	Evaluate clinical implementation and perceptions of MCEd test <i>Enrollment: complete, published</i>	<i>The Lancet 2023</i>
3	SYMPLIFY (n=6,242)	Assess MCEd test in individuals with signs/symptoms of cancer <i>Enrollment: complete, published</i>	<i>Lancet Oncology 2023</i>
4	NHS-GALLERI (n≈142,321)	Assess clinical utility of MCEd for population screening in the UK <i>Enrollment: complete</i>	
5	STRIVE (n=99,481)	Evaluate MCEd test performance in women to detect invasive cancers^a <i>Enrollment: complete</i>	
6	SUMMIT (n=13,035)	Clinical validation in individuals at high risk of lung cancer <i>Enrollment: complete</i>	
7	REFLECTION (n≈17,000)	Assess experience/clinical outcomes in real-world setting <i>Enrollment: ongoing</i>	
8	PATHFINDER 2 (n≈35,000)	Evaluate MCEd test performance in eligible screening population <i>Enrollment: completed</i>	
9	REACH (n≈50,000)	Understand health equity impact of Galleri in a Medicare population <i>Enrollment: ongoing</i>	



Cancer types detected in Substudy 3 of CCGA


45 cancer types lack recommended screening

Breast	Lung	Cervical	Colorectal Prostate
<p>Adrenal Cortical Carcinoma</p> <p>Ampulla of Vater</p> <p>Anus</p> <p>Appendix, Carcinoma</p> <p>Bile Ducts, Distal</p> <p>Bile Ducts, Intrahepatic</p> <p>Bile Ducts, Perihilar</p> <p>Bladder, Urinary</p> <p>Bone</p> <p>Esophagus and Esophagogastric Junction</p> <p>Gallbladder</p> <p>Gastrointestinal Stromal Tumor</p> <p>Gestational Trophoblastic Neoplasms</p> <p>Kidney</p> <p>Larynx</p> <p>Leukemia</p>	<p>Liver</p> <p>Lymphoma (Hodgkin and Non-Hodgkin)</p> <p>Melanoma of the Skin</p> <p>Merkel Cell Carcinoma</p> <p>Mesothelioma, Malignant Pleural</p> <p>Nasal Cavity and Paranasal Sinuses</p> <p>Nasopharynx</p> <p>Neuroendocrine Tumors of the Appendix</p> <p>Neuroendocrine Tumors of the Colon and Rectum</p> <p>Neuroendocrine Tumors of the Pancreas</p> <p>Oral Cavity</p> <p>Oropharynx (HPV-Mediated, p16+)</p> <p>Oropharynx (p16-) and Hypopharynx</p> <p>Ovary, Fallopian Tube and Primary Peritoneum</p> <p>Pancreas, exocrine</p> <p>Penis</p>	<p>Plasma Cell Myeloma and Plasma Cell Disorders</p> <p>Small Intestine</p> <p>Soft Tissue Sarcoma of the Abdomen and Thoracic Visceral Organs</p> <p>Soft Tissue Sarcoma of the Head and Neck</p> <p>Soft Tissue Sarcoma of the Retroperitoneum</p> <p>Soft Tissue Sarcoma of the Trunk and Extremities</p> <p>Soft Tissue Sarcoma Unusual Histologies and Sites</p> <p>Stomach</p> <p>Testis</p> <p>Ureter, Renal Pelvis</p> <p>Uterus, Carcinoma and Carcinosarcoma</p> <p>Uterus, Sarcoma</p> <p>Vagina</p> <p>Vulva</p>	


Galleri does not detect all cancers and all cancers cannot be detected in the blood. False positive and false negative results do occur.

PATHFINDER


Cancer signal was detected in 1.4% (92/6621 participants)



MCED Detected Cancers



Cancer Signal Origin May Aid in the Time to Diagnostic Resolution



PATHFINDER true positive participants had a faster time to diagnosis than previous retrospective analysis of cancer patients

Consistent Results Across Studies

Clinical Validation Study (CCGA3)

0.5%
False
positive rate

44%
Positive predictive value

89%
Localization accuracy

Confirmatory Intended Use Population Study (PATHFINDER)*

0.5%
False positive rate

43%
Positive predictive value


88%
Localization accuracy**

Refined test used commercially; **1st or 2nd location prediction

Klein EA, et al. Ann Oncol. 32:1167, 2021
Schrag et al. Lancet 402:1251, 2023


Modeled Stage Shift

Expected Stage-Shift in Incident Cancers by Adding the Galleri test to SoC, based on Modeled Data



Based on modeled data of GRAIL's MCED test in elevated risk population age 50–79 years.

Stage IV Cancers Account for an Outsized Proportion of Cancer Deaths




Source: SEER, adults aged 50-79, 2006-2015

Clarke et al.: Cancer Epidemiol Biomarkers Prev. 2020 29:895-902

Potential for Earlier Detection to Save Lives


Modeled Data from SEER and CCGA




**26% Reduction in
Cancer Mortality**

Methylated DNA Biomarkers and Incident Cancer in the American Cancer Society (ACS) Cancer Prevention Study-3 (CPS-3) Cohort

Supports 1-year screening interval




 **323 days**
average time between blood draw and diagnosis among those detected

Cohort sampling design
1:1 cancer:noncancer (N = 1,425 each)
Total eligible N = 245,171


High Percentage of Stage I Diagnoses With Repeat Testing

Repeat testing may improve early detection of multiple cancer types, including those currently without USPSTF-recommended screening^{a,b}





Projected Increased Relief, Happiness, and Confidence After an NCSD MCED Test Result: General Population




Relief was the most important emotional or psychological impact

Participants anticipated feeling more relief, happiness, confidence, and motivation

Participants anticipated feeling less concern/worry, stress, anxiety, fear, depression, and frustration



Impacts of NCSD MCED Test Result Include Positive Changes in Health-Related Behaviors: General Population



Majority (73%) of participants would maintain or improve their health-related behaviors

Most participants were very-to-extremely certain of increasing health-related behaviors, including adherence to recommended cancer screening and eating healthy

Despite this

USPSTF Recommendations for Cancer Screening

Cancer	Grade	Population	Modality/ Recommendation
Cervical	A	Women aged 21 to 65	Regular screening (3–5 years) using cervical cytology and/or HPV tests
Colorectal	A B	Adults aged 50 to 75 Adults aged 45-49	Regular annual screening, multiple effective methods available
Breast	B C	Women aged 50 to 74 Women aged 40 to 49	Biennial screening mammography
Lung	B	Adults aged 55–80, with history of smoking	Annual low-dose computed tomography (LDCT) screening
Prostate	C	Men aged 55 to 69	Periodic PSA screening on case-by-case basis




Mortality



> 600,000 people die of cancer every year In the US

To achieve this

Adding MCED has the potential...



26% Reduction in Cancer Mortality

Characteristics



Female



67



No specific risk factors noted beyond age-associated risk and maternal history of multiple myeloma

Stage IA Ovarian Clear Cell Carcinoma

Pathological Diagnosis Confirmed at 28 Days Following MCED Test Results

Pelvic Ultrasound and CT Scan

- 8.5x8.7cm mass of suspected uterine origin.
- Subsequent CT scan confirmed vascular pelvic mass arising from left ovary.

Surgery

- Total abdominal hysterectomy, bilateral salpingo-oophorectomy, bilateral pelvic lymph node dissection, and peritoneal biopsy and washing.

Treatment

- Treatment with 6 cycles of adjuvant chemotherapy (carboplatin and taxol) initiated.

Blood Draw




MCED test report returned

Cancer Signal Detected

CSO Prediction

CSO1=Uterus
CSO2=Ovary



Day 1

Day 10

Day 28

Day 64

Diagnosis

- Surgical pathology showed a stage pT1AN0M0 ovarian clear cell carcinoma.


Day 28

Current Status

- A 6-month follow up CT scan was negative, and individual remains asymptomatic with **no clinical evidence of disease after 21 months of follow up.**

NED
21
Months

STATUS QUO BIAS



WOULD YOU
LIKE A
SHOVEL?



... BUT I'M
ALREADY USING
THIS SPOON!

