

Multi-Cancer Detection Test to Aid Clinically Presenting Head and Neck Cancer Diagnosis

American Academy of Otolaryngology-Head and Neck Surgery (AAO-HNSF) 2021 Annual Meeting & OTO Experience, October 3-6, 2021

Alan H. Bryce,¹ Minetta C. Liu,² Michael V. Seiden,³ David D. Thiel,⁴ Donald Richards,⁵ Carlos Becerra,⁵ Kathryn N. Kurtzman,⁶ Xiaoji Chen,⁶ Tony Wu,⁶ Quan Zhang,⁶ Jingjing Gao,⁶ Nan Zhang,⁶ Earl Hubbell,⁶ Arash Jamshidi,⁶ Eric T. Fung,⁶ Eric A. Klein⁷
¹Mayo Clinic, Phoenix, AZ; ²Mayo Clinic, Rochester, MN; ³US Oncology Research, The Woodlands, TX; ⁴Mayo Clinic Florida, Jacksonville, FL; ⁵US Oncology Research, Tyler, TX; ⁶GRAIL, LLC, a subsidiary of Illumina, Inc., Menlo Park, CA; ⁷Cleveland Clinic, Cleveland, OH

INTRODUCTION

- Delayed diagnostic resolution in individuals seeking care for signs and symptoms suspicious for cancer can affect outcome, patient experience, and the efficiency of care delivery.¹⁻⁶
- The Circulating Cell-free Genome Atlas (CCGA; NCT02889978) study was designed to develop and validate a multi-cancer detection test.
- This blood-based multi-cancer detection test uses a cell-free DNA (cfDNA)-based targeted methylation assay to detect cancer signals across multiple cancers, including head and neck cancer, and predict cancer signal origin (ie, the tissue where the primary tumor arose).⁷

OBJECTIVE

- To evaluate performance of the multi-cancer detection test in a subgroup of clinically presenting cancer participants, including a smaller subgroup of head and neck cancer, and non-cancer participants.

RESULTS

Study Population

- A total of 518 clinically presenting cancer participants, 25 of which had head and neck cancer were included in the analysis. A subgroup of 327 clinically presenting cancer participants with single known primaries and a trained, predicted CSO (not including a CSO prediction of “other”) were included in the accuracy of CSO prediction analysis.
- In addition, 398 non-cancer participants with confirmed non-cancer status at Year 1 follow-up were analyzed. Also, 16 participants who had significant non-malignant medical conditions at enrollment were analyzed.

Specificity

- The multi-cancer detection test correctly identified 396 of 398 non-cancer participants (non-cancer status confirmed at Year 1 follow-up) as not having cancer (**Table 1**), corresponding to a false positive rate of <1%.
- High specificity was also observed in non-cancer participants with non-malignant confounding conditions.

Sensitivity

- Sensitivity of the multi-cancer detection test in participants with clinically presenting cancers and the subgroup of participants with head and neck cancers was comparable (**Table 2**).
- Sensitivity was generally higher for head and neck cancers at later clinical cancer stages, as expected.

- It should be noted that the wide confidence intervals may be attributed to small sample size.

Cancer Signal Origin Prediction

- Accuracy of CSO prediction for participants with head and neck cancers was high, and similar to that of all clinically presenting cancer participants (**Table 3**).

Table 1. Specificity

Group	Cancer Signal Not Detected/Total	% (95% CI)
Non-cancers with non-cancer status confirmed at Year 1 follow-up	396/398	99.5% (98.2-99.9%)
Non-cancers with non-malignant confounding conditions at enrollment	15/16	93.8% (71.7-99.7%)

CI, confidence interval.

Table 2. Sensitivity

Group	Cancer Signal Detected/Total	% (95% CI)
All Clinically Presenting Cancer	344/518	66.4% (62.2-70.3%)
Clinically Presenting Head and Neck Cancer	18/25	72.0% (52.4-85.7%)
Stage I	2/3	66.7% (20.8-98.3%)
Stage II	1/4	25.0% (1.3-69.9%)
Stage III	5/6	83.3% (43.6-99.1%)
Stage IV	10/12	83.3% (55.2-95.3%)

CI, confidence interval.

Table 3. Cancer Signal Origin Prediction Accuracy

Group	Correct Signal Origin Detected/Total	% (95% CI)
All Clinically Presenting Cancer ^a	300/327	91.7% (88.3-94.3%)
Clinically Presenting Head and Neck Cancer	16/18	88.9% (67.2-96.9%)

^aExcludes participants with multiple or unknown primaries or other cancers. CI, confidence interval.

CONCLUSIONS

- The multi-cancer detection test performed with high specificity while detecting cancer signals and predicting cancer signal origin in individuals with symptomatic cancer, including in those with head and neck cancer.
- These findings support the potential use of this multi-cancer detection test for diagnostic resolution of symptomatic cancers.

References

- Jensen H, et al. *BMC Cancer*. 2014;14:636.
- Lyratzopoulos G, et al. *Br J Cancer*. 2015;112(suppl 1):S84-S91.
- Neal RD, et al. *Br J Cancer*. 2015;112(suppl 1):S92-S107.
- Mendonca SC, et al. *Eur J Cancer Care (Engl)*. 2016;25(3):478-490.
- Chapman D, et al. *Br J Cancer*. 2020;123(5):722-729.
- World Health Organization. Guide to cancer early diagnosis. World Health Organization; 2017. <https://apps.who.int/iris/handle/10665/254500>
- Liu MC, et al. *Ann Oncol*. 2020;31(6):745-759.

Disclosures

AHB has consulted for Pfizer, Merck, Bayer, and Astellas Pharmaceuticals. MCL is an uncompensated advisory board member of GRAIL, LLC, a subsidiary of Illumina, Inc. MVS is a full-time employee of and has stock in McKesson Corporation and serves as a clinical advisor for GRAIL, LLC, a subsidiary of Illumina, Inc. DDT, DR, and CB have no conflicts of interest to declare. KNK, TW, QZ, AJ, ETF, KNK, TW, QZ, AJ, ETF, XC, JG, JG, NZ, EH are full-time employees of GRAIL, LLC, a subsidiary of Illumina, Inc., with equity in Illumina, Inc. EAK is a consultant of GRAIL, LLC, a subsidiary of Illumina, Inc.

Acknowledgements

Funded by GRAIL, LLC, a subsidiary of Illumina, Inc. Writing and editorial assistance provided by Prescott Medical Communications Group (Chicago, IL).

METHODS

Multi-Cancer Detection Test

- This test detects cancer signal and predicts cancer signal origin through machine-learning analysis of methylation patterns in cfDNA circulating in the blood.⁷
- Targeted methylation sequencing of 10⁶ CpG sites in cfDNA isolated from a blood draw is performed.
- A machine-learning algorithm differentiates abnormal tumor cfDNA methylation patterns from normal cfDNA methylation patterns.
- When a cancer signal is detected, the machine-learning algorithm discerns cancer type-specific differences in methylation to predict cancer signal origin.

Study Design and Participants

- The CCGA study is a prospective, longitudinal, multicenter, case-control study with the objective to develop and validate a multi-cancer detection test.
- In the second CCGA substudy, the multi-cancer detection test performance was assessed in a training set and an independent validation set.⁷
- Test performance is reported from the validation set from the second CCGA substudy of participants.
- Specificity was assessed in non-cancer participants who had confirmed non-cancer status at Year 1 follow-up and in a subgroup of non-cancer participants with significant non-malignant medical conditions at enrollment (eg, oral/pharyngeal/nasal lesion).
- Sensitivity and cancer signal origin prediction accuracy were assessed in a subgroup of participants with clinically presenting cancers and a smaller subgroup of participants with head and neck cancer from the second CCGA substudy.

Measures

- **Specificity:** Estimated as the proportion of non-cancer participants with a “cancer signal not-detected” result.
- **Sensitivity:** Estimated as the proportion of cancer participants with a “cancer signal detected” result.
- **Sensitivity by clinical stage:** Estimated as the proportion of participants with a “cancer signal detected” result among cancer participants with a clinical stage.
- **Accuracy of cancer signal origin (CSO) prediction:** Estimated as the proportion of participants with a true positive cancer signal detected result and with a correctly predicted trained cancer signal origin.